


## Produced for the Optometry Blog by CLEARVIEW Training

### For best results:

1. In your printer menu make sure that the option of 'Print Scaling' is set to 'None.'
2. For maximum contrast, use an Inkjet printer and set to 'Best Quality' in your Properties in the Printer menu. The contrast will not be high enough if you use a LaserJet printer
3. When printed the height of the letters should measure 87.3 mm
4. Laminate the letters with a matt laminate pouch to avoid veiling glare.



6/60 (87.3mm)

**N**

6/60 (87.3mm)



6/60 (87.3mm)